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Cat: SB-WB011; **Size:** 50 ml+50 ml; 5*(50 ml+50 ml)

Product Description

The LumiBest ECL Substrate is an excellent luminescent substrate developed for HRP systems, which can be used for HRP-labeled Western Blot and ELISA experiments. X-ray film or luminometer can be used to detect the light intensity or image. The light intensity of the labeled antigen will be qualitatively or quantitatively detected. LumiBest has revolutionized the traditional ECL substrate system, resulting in excellent detection limits, duration of luminescence and reagent stability. Overall performance evaluation parameters would be enhanced more than one billion times. LumiBest exceed the world-renowned fg-level sensitivity systems in actual performance.

LumiBest has the following features

Sensitive: high S/N ratio, low background, femtomole antigens detectable.

Lasting: Provides plenty of time to get the depth of the desired band.

Convenient: Bands immediately apparent without waiting.

Stable: Long-term storage does not reduce the sensitivity. **Compatibility:**

Fit the operating habits from different laboratories. Strong compatibility is easy for novices to get started as well as easy for experts to replace the old system without optimization.

Product Character:

Liquid A: light yellow to colorless liquid; Liquid B: colorless liquid

Materials needed

1. Secondary antibody incubated Western blot membrane.
2. TBST or PBST buffer.
3. Chemiluminescence detection equipment or X-ray film.
4. If using X-ray film you need cassette and plastic film.

How to use

1. Western Blot experiments were performed as usual and the secondary antibody incubation was completed. The membrane was washed three times for 5 min each with TBST or PBST buffer.
2. Take a certain volume of liquid A and liquid B, the 1:1 mixing reserve (For example, a 6CM * 9CM PVDF membrane needs 400ul + 400ul liquid A + B, skilled experimenter can also add less volume).
3. Remove the TBST or PBST buffer in the container. The mixed liquid A + B were evenly leached on the membrane.

Troubleshooting

Problem	Cause	Solution
Yellow band	Secondary antibody concentration is too high	Reduce the secondary antibody concentration
Light-emitting time getting shorter		
Black band		
Weak signal or no signal	Low target protein content in Sample	Add more antigen
	No target protein in the sample	Check up sample preparation
	Incomplete transferring	Extend the transfer time
	Second antibody degeneration	Make new secondary antibody solution
Background is too dark or appear non-specificity bands	Poor primary antibody specificity	Adjusted primary antibody concentration
	Incomplete blocking	Extend the blocking time
	Blocking solution degeneration	Make new blocking solution
	Overexposure	Reduces exposure time